





Annexe 9: Skills Checklist: Weighing the Baby


Scale Type 1:

| Picture/Illustration | Skills Checklist | For Practice | | | | |
|---|---|--------------|---|---|---|---|
| | | 1 | 2 | 3 | 4 | 5 |
|  | 1. Place the sling on scale | | | | | |
| | 2. Hold scale by top bar off the floor, keeping the adjustment knob at eye level. | | | | | |
| | 3. Turn the screw until its top fully covers the red and 'O' is visible. | | | | | |
|  | 4. Remove sling on hook and place it on a clean cloth on the ground. | | | | | |
| | 5. Place body with minimum clothes on, in sling and replace the sling on hook. | | | | | |
|  | 6. Holding top bar carefully, as you stand up, lift the scale and sling with baby off the ground, until the knob is at eye level. | | | | | |
|  | 7. Read the weight. | | | | | |
| | 8. Gently put the sling with baby in it, on the ground and unhook the sling. | | | | | |
| | 9. Remove the baby from the sling and hand it over to its mother. | | | | | |
| | 10. Record the weight. | | | | | |

Scale Type 2:



Salter Scale



Weighing machine (for up to 10 kg)